

DPA ANAESTHESIA SYSTEM SPECIFICATION



DPA

- The Diamedica Portable Anaesthesia system is:
- 47cm width, 33cm height, 17cm depth and weighs 9.5 kg.
- Designed for field use it will accept any oxygen source available.
- Ventilation is carried out manually using an in-line self-inflating bag.



VAPORISER

- Low resistance vaporiser suitable for drawover anaesthesia
- Suitable for continuous flow anaesthesia
- Stainless steel construction
- Agents halothane and isoflurane
- Scale 0 to 5%
- Capacity 150ml
- Weight 2.6Kg empty
- Maintenance none unless thymol from halothane builds up. Can be cleared using fresh halothane in a few seconds

